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**Steriles Produkt zum Einmalgebrauch**  
**Monopolare HF-Elektrode Art.-Nr.:** 8M910381 (205-0s), 8M910382 (210-3s), 8M910383 (214-2s), 8M910376 (300-0s), 8M910377 (303-2s), 8M910378 (305-4s), 8M910380 (305-6s)

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**1 Scope**

**Sterile Product for single use**

**Monopolar HF Electrode Art. No.:** 8M910381 (205-0s), 8M910382 (210-3s), 8M910383 (214-2s), 8M910376 (300-0s), 8M910377 (303-2s), 8M910378 (305-4s), 8M910380 (305-6s)

**Maximum rated voltage of accessory (U<sub>max</sub>):**

Description	U <sub>max</sub>
Monopolar HF Electrodes	5,0 kVp

**See also labeling or catalogue.**

In any combination with another electrosurgical accessory, the maximum rated voltage of the combination corresponds to the lowest rated voltage of the accessories used. (See also section "Intended Use")  
See these instructions for use, the label, or the current product catalog for the maximum rated voltage of the product.

Prior to usage, read the entire IFU of this product and of any used accessory as well as HF-generator and HF-neutral electrode (monopolar application).  
All requirements, safety notices and warnings included in the respective IFUs have to be followed strictly.  
This medium is not intended for users in the USA.

**2 Intended Use**

Product intended for use in open or rather endoscopic surgery.

**For cutting and coagulating of biological tissue.**

Only to be used by skilled medical professional who have been introduced to the product, for example by this IFU.

Connect product to the monopolar output of the HF-generator by using the integrated cable of the device.

The product has to be inserted in a suitable HF-Cable or HF-Handle, activation is done by the buttons of the handle or with the footswitch.

It is strictly prohibited to manipulate the electrode tip.

The product is supplied sterile and is for single use only.

It is strictly prohibited to reuse or reprocess the product.

**Indication:**

Product intended for use in open or rather endoscopic surgery.

**Contraindication:**

Product is not intended for direct contact with the heart or the central circulatory or central nervous system.

Application of high-frequency current may interfere with cardiac pacemakers and in vivo heart defibrillators, so affected patients must consult a cardiologist prior to the intervention.

It is not allowed the patient contacts any grounded metal objects during the intervention.

**Connection and Activation:**

Prior to application of monopolar products, it has to be safeguarded that an HF-neutral electrode has been applied correctly to the patient and that this device is connected properly to the respective HF-generator.

Activation is done using the buttons of the HF-handle or with the foot switch of the HF-generator.

**Combination / Compatibility:**

Prior to use, compatibility of the product with the intended HF-generator has to be verified.

In case of uncertainties, contact the manufacturer of this product or the manufacturer of the used HF-generator.

The frequency of the used HF-generator shall not exceed the maximum frequency of 4 MHz and not exceed the maximum rated voltage of the accessory (see section "Scope").

Follow instructions, safety notices and warnings included in the User Manuals or IFU of the used accessories and used HF-generator.

**It is recommended to use a smoke evacuation system.**

**IMPORTANT:**

Handle with utmost care.

This does not only apply for the duration of the surgery but also for the complete duration of storage and transport as well as during the process of connecting the product with the respective HF-handle, HF-cable and/or HF-generator.

This applies especially for the thin tips and other sensitive areas, for example the insulation.

Follow instructions and warnings in the IFU of HF-generator.

**Side effects:**

Electrosurgical interventions may lead to undesired damages or undesired burns of surrounding tissue.

Tissue damage due to extensive forces during butt dissection is possible.

It is the sole responsibility of the physician to include this aspect into their professional decision regarding the selection of treatment process.

All electrosurgical products potentially may lead to muscle stimulation during the application.

Improper use immediately will result in loss of warranty.

Liability for any damages incurred will not be accepted.

**3 Safety Notices - WARNING**

See these instructions for use, the label, or the current product catalog for the maximum rated voltage of the product.

Make sure the User Manual of the HF-Generator has been read carefully before using the product.

Prior to usage, a visual inspection and function test has to be done (see section "Visual Inspection and Function Test")

In case damages, deformation or similar is detected on the product, it is not allowed to use the device.

It has to be replaced by a new one.

At least one (1) sterile backup product has to be available for the application.

It is the responsibility of the user to determine the appropriate product size and product type according to their professional judgement and based on the patient's specific indication, preferred surgical technique and history, etc.

Prior to application, verify that the product is inserted firmly in the HF-handle/HF-cable.

This must be done carefully, in order to avoid damages on the product and/or injuries of patient, surgical personnel or third party.

Excessive force can damage the product.

Therefore, the product has to be observed during the complete application.

Verify damage free sterile packaging and expiration date

Do not use a product that is included in a damaged sterile packaging or that is beyond the expiration date.

Maintenance of the sterile condition after sterilization process is the sole responsibility of the user.

Prior to usage a visual inspection has to be done for pressure points of damages on the product.

**Exclusion:**

Do not activate the HF-generator during connection of the product.

As long as the product is activated, do not insert the tip in the surgical site and do not remove it when is activated.

Do not activate the product as long as it is in contact with metal objects and/or optics.

During an electrosurgical intervention the patient must not come into contact with grounded metal objects such as surgical desk frames, instrument trays etc.

"Do not use in the product in the immediate vicinity of inflammable or combustible materials or substances. Active electrodes must not come into contact with inflammable or combustible substances such as disinfectants, gases etc. as otherwise there is a danger of fire or explosion due to heat or sparking."

Do not activate the product uninterruptedly over a longer period.

**(Refer to section "General Safety Notices for HF-Technology")**

After switching off the electrosurgical current, the product tip, still may be hot, so that this could lead to undesirable burns.

While the product is activated, do not insert or withdraw the tip from the surgical site.

Inadvertent activation of the product or moving the tip outside of the field-of-vision can lead to unintentional injury of tissue.

Do not use metal cannula, as these may damage the product insulation (capacitive coupling). This could lead to undesirable burns.

It is recommended to use cannula systems made of plastics.

It is the responsibility of the user to apply low performance settings on the HF-Generator in order to achieve the envisaged effect for the respective intervention.

During the surgery, mechanical forces may lead to deformation or signs of wear of the product.

If, abrasion of the product is detected during the application, the product has to be exchanged by a new one.

Abrasion, adhesion of tissue, discoloration, sooting etc. do not represent a reason for complaint and do not permit claiming the manufacturer's warranty.

**4 General Safety Notices for HF-Technology (Excerpt)**

In addition to the acknowledged benefits of HF-Technology, the application includes several risks that have to be attended:

Improper use and non-observance of the IFU can lead to unintentional burns of the patient, user or third party.

Continued further education of the surgical personnel is recommended.

**a) Environment**

Pay attention that no flammable substances (anesthetics, oxidizing gases, endogenous gases etc.) are present in immediate vicinity as otherwise there is danger of explosion.

Only use non-flammable cleaning or disinfection agents.

All oxygen connections must be tight and leak-proof during the procedure.

**b) Patient Positioning and Preparation**

Ensure proper patient positioning, i.e., only use insulating surgical drapes that are dry, absorbent and liquid-tight.

Isolate conductive surfaces and contact points towards the patient.

Dry tissues of pulp are necessary to use for skin and breast folds as well as between the extremities.

Prior to application, remove any liquids that potentially accumulated in body cavities

Only use non-flammable disinfectants.

Do not use alcohol-based tinctures.

Only use non-conductive rinsing fluid, if medicinally possible.

Attend the requirements on irrigation fluid for monopolar and bipolar products.

Remove all body jewelry from the patient.

Putting a band aid over the body jewelry is not sufficient!

**c) Connections**

Prior to application, make sure the used HF-handle/HF-cable is connected properly to the HF-generator and make sure the correct power and performance setting is adjusted and displayed.

It is the responsibility of the user to apply low performance settings on the HF-Generator in order to achieve the envisaged effect for the respective intervention.

**d) HF-Neutral Electrode for Monopolar Application**

In case of monopolar application, select a HF-neutral electrode suitable for the patient, apply it correctly and connect it properly with the respective HF-generator.

Follow all instructions for proper application of the HF-neutral electrode, incl. patient protection and patient monitoring, monitoring of the HF-neutral electrode and all further provisions, safety notices and warnings included in the IFU of the HF-neutral electrode.

**e) Patient Reactions**

All electrosurgical devices potentially can cause muscle stimulation during the application.

The design of this product minimizes the risk of this undesirable effect.

Nevertheless, muscle stimulation can lead to an unexpected movement of the patient in the surgical field.

**f) Handling HF-Accessory**

Make sure the accessory is compatible.

Do not touch the instrument tip during the complete application.

As long as the product is not applied, place it on a dry, clean and non-conductive and well-visible surface, that is not in contact with the patient.

Never store product on the patient.

Unintended activation of the product can lead to unintentional burns or other injuries of patient, user or third party.

Keep on longer breaks between activation phases.

Only adjust low power settings.

**g) Completeness of the System**

At the end of the surgery, confirm completeness of the system.

**5 Visual Inspection and Function Test**

Prior to usage, check the sterile packaging for damage and the expiry date on the sterile packaging.

A product whose sterile packaging is damaged and/or whose expiry date has passed must not be used.

It must be replaced with a new one.

Never place the product on the patient or in the immediate vicinity of the patient.

Do not use damaged product.

Prior to usage, check the entire product, especially the insulation and product tip, for pressure marks and damage.

A product with damage, pressure marks or questionable condition must not be used and must be replaced with a new one.

During and after use, tissue build-up or sooting may occur on the distal end of the active electrode.

Such adhesions or sooting do not constitute grounds for complaint; the product must be replaced with a new one.

Due to a longer period of use, mechanical forces or similar, deformation or wear of the insulation material may occur.

This is also not a reason for complaint, the product must be replaced with a new one.

A blockage of the suction channel (if present) is not a reason for complaint.

A product whose suction channel (if any) is clogged must be replaced with a new one.

An electrical continuity test must be carried out before use.

If the electrical continuity test is not passed, the product must not be used and must be replaced with a new one.

This is also not a reason for complaint, the product must be replaced with a new one.

A blockage of the suction channel (if present) is not a reason for complaint.

A product whose suction channel (if any) is clogged must be replaced with a new one.

An electrical continuity test must be carried out before use.

If the electrical continuity test is not passed, the product must not be used and must be replaced with a new one.

**6 Exclusion of Repair and Modification**

Unauthorized repairs and modification are strictly prohibited and lead to an immediate loss of manufacturer's warranty.

It is not allowed to bend product or to manipulate the product in any other manner.

Especially product having a hook as active part is never allowed to be bent.

This could lead to severe injuries of patient, user or third party.

The product is delivered sterile and for single use only.

It is not allowed to reuse or reprocess the product.

See also labelling.

**7 Packaging and Labelling**

Only use product in closed original packaging with an undamaged labelling.

In case of uncertainties, contact the manufacturer.

Only use product of an undamaged sterile packaging and valid expiration date.

Do not use product, that exhibits a damaged sterile packaging or that is beyond the expiration date.

See also labelling.

**8 Storage**

Store product in a clean and dry environment protected from direct sunlight.

Store product in the closed original packaging.

Maintenance of the sterile condition after sterilization process is the sole responsibility of the user.

Do not use the product any longer if the expiration date is expired or that has a damaged sterile packaging.

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info@dr-medical.ch | www.dr-medical.ch

**11 Warranty**

The products satisfy the highest quality standards.

Liability and warranty are excluded for all products that have been modified in any way, not applied according to their intended use, or that have been handled or applied improperly, or in case of any other deviation from the instructions set forth in this IFU.

Furthermore, the manufacturer denies any liability for any accidental, intentional damage or for a damage or loss arising out of handling or application of the product.

Additionally, all liability and warranty are omitted in case our product was repaired by a company that has not been authorized by us.

Unauthorized repairs are strictly prohibited.

**12 Return**

Returned products only are accepted, if they are marked as "hygienically safe" or "not contaminated" and have been packed safely for the transport.

Use our return sheet for any return of products.

**13 Disposal**

Disposal of the product, their packaging material as well as any accessory has to be done according to the applicable country specific requirements, regulations and laws.

Additionally, the applicable requirements of the respective clinical place(s) in regards to disposal of medical devices have to be followed as well.

**14 Regulatory Remark**

Due regulatory reasons we would like to inform users and patient, that serious incidents in connection with our medical device have to reported to the manufacturer and competent authority of the EU Member State where user and/or patient are resident.

**15 About these Instructions for Use**

Throughout the period of usage these IFU must be kept freely accessible for the user.

For a current revision of these IFU, please contact medimex GmbH or visit the medimex Homepage for a download of this document:

www.medimex.de/downloads

Changes reserved.

This medium is not intended for users in the USA.

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REGER Medizintechnik GmbH, Gewerbestrasse 10, 78667
DR MEDICAL AG, Weissensteinstr. 26, Ziegelmatteareal Bau 6, 4503 Solothurn