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1 Scope

Non-sterile, single-use product
Monopolar HF Electrodes, Article No.: 205-5SU, 210-3SU, 214-2SU
Maximum rated voltage of accessory (U_{max}):

Article No.	U _{max}
205-5SU, 210-3SU, 214-2SU	4,3 kVp

See also labeling or catalogue.

In any combination with another electrosurgical accessory, the maximum rated voltage of the combination corresponds to the lowest rated voltage of the accessories used. (See also section "Intended Use")
See this Instructions for Use, the label, or the current product catalog for the maximum rated voltage of the product.
In case of uncertainties, contact the manufacturer.
Prior to usage, read the entire IFU of this product and of any used accessory as well as HF-generator.
All requirements, safety notices and warnings included in the respective IFUs have to be followed strictly.
This medium is not intended for users in the USA.

2 Intended use

This product is intended for open or endoscopic surgery and is used for cutting and coagulating biological tissue.

Only for use by skilled medical professionals.

Indication:
Product intended for use in open or rather endoscopic surgery.

Contraindication:

Product is not intended for direct contact with the heart or the central circulatory or central nervous system.

Application of high-frequency current may interfere with cardiac pacemakers and in vivo heart defibrillators, so affected patients must consult a cardiologist prior to the intervention.

Connection and Activation:

Prior to application of monopolar products, it has to be safeguarded that an HF-neutral electrode has been applied correctly to the patient and that this device is connected properly to the respective HF-generator.

The product is connected by a suitable HF-cable or a suitable HF-handle to the respective exit of an HF-generator.

Activation is done using the buttons of the HF-handle or with the foot switch of the HF-generator.

Combination/Compatibility:

Prior to use, compatibility of the product with the intended HF-handle/HF-cable and HF-generator has to be verified.

The frequency of the used HF-generator shall not exceed the maximum frequency of 4 MHz and not exceed the maximum rated voltage of the accessory (see section "Scope").

In case of uncertainties, contact the manufacturer of this product or the manufacturer of the used accessory or HF-generator.

Follow the instructions, safety notices and warnings included in the User Manuals or IFU of the used accessories and used HF-generator.

It is recommended to use a smoke evacuation system.

IMPORTANT:

Handle with utmost care.

This does not only apply for the duration of the surgery but also for the complete duration of storage, processing and transport as well as during the process of connecting the product with the HF-accessory and HF-generator.

This applies especially for the thin components and other sensitive areas, e.g. the insulation.

In case of use, the product has to be cleaned, disinfected and sterilized (see section "Processing: Cleaning, Disinfection and Sterilization").

Improper use immediately will result in loss of warranty.

Liability for any damages incurred will not be accepted.

3 Safety notices – WARNING!

It is the responsibility of the user to determine the appropriate product size according to their professional judgement and based on the patient's specific indication, preferred surgical technique and history, etc.

See these instructions for use, the label, or the current product catalog for the maximum rated voltage of the product.

Prior to application, the product has to be cleaned, disinfected and sterilized according to a validation procedure (EN ISO 17665).

Prior to use a visual inspection and function test has to be done (see section "Visual inspection and function test").

In case of damages, deformation or similar is detected on the product, it is not allowed to use the device.

It has to be replaced by a new one.

At least one (1) cleaned, disinfected and sterilized backup product has to be available.

It is the responsibility of the user to determine the appropriate product size and product type according to their professional judgement and based on the patient's specific indication, preferred surgical technique and history, etc.

Prior to application, verify that the product is inserted firmly in the HF-handle/HF-cable.

This must be done carefully, in order to avoid damages on the product and/or injuries of patient, surgical personnel or third party.

Excessive force can damage the product.

Therefore, the product has to be observed during the complete application.

Exclusion:

Do not activate the product as long as it is in contact with metal objects and/or optics.

During an electrosurgical intervention the patient must not come into contact with grounded metal objects such as surgical desk frames, instrument trays etc.

Pay attention that no flammable substances are present in immediate vicinity as otherwise there is danger of explosion.

Do not use in the product in the immediate vicinity of inflammable or combustible materials or substances. Active electrodes must not come into contact with inflammable or combustible substances such as disinfectants, gases etc. as otherwise there is a danger of fire or explosion due to heat or sparking.

Do not activate the product uninterruptedly over a longer period.

(Refer to section "General Safety Notices for HF-Technology")

After switching off the electrosurgical current, the product tip may still be hot and can lead to unintentional burns.

Do not activate the product uninterruptedly during a longer period.

It is the responsibility of the user to apply low performance settings on the HF-Generator in order to achieve the envisaged effect for the respective intervention.

During the surgery, mechanical forces may lead to deformation or signs of wear of the product. If, abrasion of the product is detected during the application, the product has to be exchanged by a new one.

Abrasion, adhesion of tissue, discoloration, sooting etc. do not represent a reason for complaint and do not permit claiming the manufacturer's warranty.

4 General Safety Notices for HF Technology (Excerpt)

In addition to the acknowledged benefits of HF-Technology, the application includes several risks that have to be attended:

Improper use and non-observance of the IFU can lead to unintentional burns of the patient, user or third party.

Continued further education of the surgical personnel is recommended.

a) Environment

Pay attention that no flammable substances (anesthetics, oxidizing gases, endogenous gases etc.) are present in immediate vicinity as otherwise there is danger of explosion.

Only use non-flammable cleaning or disinfection agents.

All oxygen connections must be tight and leak-proof during the procedure.

b) Patient Positioning and Preparation

Ensure proper patient positioning, i.e., only use insulating surgical drapes that are dry, absorbent and liquid-tight.

Isolate conductive surfaces and contact points towards the patient.

Dry tissues of pulp are necessary to use for skin and breast folds as well as between the extremities.

Prior to application, remove any liquids that potentially accumulated in body cavities.

Only use non-flammable disinfectants.

Do not use alcohol-based tinctures.

Only use non-conductive irrigation fluid, if medicinally possible.

Attend the requirements on irrigation fluid for monopolar and bipolar products.

Prior to application, remove any body jewelry from the patient.

Putting a band aid over the body jewelry is not sufficient!

c) Connections

Prior to application, make sure the product is connected properly to the HF-generator and make sure the correct power and performance setting is adjusted and displayed.

It is the responsibility of the user to apply low performance settings on the HF-Generator in order to achieve the envisaged effect for the respective intervention.

d) HF-Neutral Electrode for Monopolar Application

In case of monopolar application, select a HF-neutral electrode suitable for the patient, apply it correctly and connect it properly with the respective HF-generator.

Follow all instructions for proper application of the HF-neutral electrode, incl. patient protection and patient monitoring, monitoring of the HF-neutral electrode and all further provisions, safety notices and warnings included in the IFU of the HF-neutral electrode.

a) Patient Reactions

All electrosurgical devices potentially can cause muscle stimulation during the application. The design of this product minimizes the risk of this undesirable effect.

Nevertheless, muscle stimulation can lead to an unexpected movement of the patient in the surgical field.

f) Handling HF-Accessory

Make sure the accessory is compatible.

Do not touch the instrument tip during the complete application.

As long as the product is not applied, place it on a dry, clean and non-conductive and well-visible surface, that is not in contact with the patient.

Never store product on the patient.
Unintended activation of the product can lead to unintentional burns or other injuries of patient, user or third party.

Never wrap cables around the patient and never lay cables over the patient.

Apply only short activation times.

Keep on longer breaks between activation phases.

Only adjust low power settings.

g) Completeness of the System

At the end of the surgery, confirm completeness of the system.

5 Processing: Cleaning, disinfection and sterilization

5.1. Maximum number of reprocessing cycles

The products mentioned in section 1 "Scope" are intended for single use only.

The reprocessing and reuse of the products is prohibited.
Visual inspection and function test, especially the condition of the insulation and product tip are decisive for whether the product can be applied.

5.2. Time requirements for cleaning and disinfection.

Avoid idle time.

5.3. Validation of Processing

The following validated processing procedure is recommended.

Equivalent deviant processes are possible.

Then, it is the sole responsibility of the user to safeguard the suitability of the actually applied procedure by suitable means (e.g. validation, routine monitoring, verification of material compatibility etc.).

Automated cleaning and disinfection always are preferred.

The following procedure was validated according to EN ISO 17665.

Additionally, applicable processing requirements specific to the respective clinical place (operator) as well as national or country specific regulations have to be followed as well.

Never use sharp objects for cleaning.

Disinfectants always have to be rinsed carefully.

5.4. Preparation for Cleaning

Remove product from their packaging.

Place it in a container provided for cleaning.

It is not necessary to disassemble the product.

5.5. Pre-Cleaning

Use tap water (potable water quality) (<40°C) and aldehyde-free, non-flxing disinfectants if applicable.

Thoroughly remove surface stain with a soft brush or synthetic fleece, as otherwise particles or dried secretions may adhere.

Ensure that areas difficult to reach are cleaned thoroughly and rinsed several times.

Cavities and lumen have to be rinsed intensively using at least 3x20 ml cold tap water (<40°C) with the aids of a rinsing adapter (e.g. from the company Medisafe), or with a syringe or with a water jet pistol (>30 Sec.).

This pre-cleaning step always has to be done prior to the manual cleaning or cleaning with the cleaning and disinfection device (CDD).

5.6. Manual Cleaning and Disinfection

Prepare an immersion bath with a suitable fluid cleaning agent.

Use a cleaning agent that is compatible to the disinfectant and suitable for immersion baths.

Follow the instructions of the manufacturer of cleaning agent and disinfectant.

Only use agents suitable for medical devices made from metal and plastics with a pH-value between 5.5 and 12.3.

Recommendation: Cleaning agent gigazyme® (Schülke & Mayr) and disinfectant Korsalex Plus.

Do not use high alkaline cleaning agents.

These will impair the lifetime of the product.

• Prepare immersion bath with cleaning agent according to the specific cleaning agent IFU.

• Prepare a separate immersion bath with disinfectant according to the specific disinfectant IFU.

• Immerse the product completely in the ultrasonic bath with cleaning agent (e.g. 0.5% gigazyme®)

• Clean product in ultrasonic bath using a sonication time of 5 Min. and a frequency of 45 kHz.

• Follow all instructions set forth in the IFU of the cleaning agent, disinfectant and ultrasonic bath.

• Ensure the product will not touch other products or parts in the ultrasonic bath.

• Prepare immersion bath with cleaning agent in the ultrasonic bath are avoided.

• Then, clean product with a soft brush under cold running town water (<40°C).

• Intensively rinse cavities and lumen with a water jet pistol (>30 Sec.) or similar for at least 1 Min.

• Afterwards, rinse product thoroughly for at least 1 min. under running tap water (>40°C) in order to remove any residues of the cleaning agent.

• Inspect product visually on remaining stain.

• In case stain is still present, repeat aforementioned cleaning steps as long as it needs until no soil is present.

• Afterwards: immerse product completely in a disinfectant bath including e.g. Korsalex Plus 3%, for at least 15 Min.

• Follow the manufacturer's data for residence time.

• Ensure the disinfectant will contact all areas on the product.

• Rinse cavities and lumen several times, that means at least 3x with 20 ml each of disinfectant bath fluid.

• Afterwards: Rinse product thoroughly for at least 1 Min. with demineralized cold water, in order to remove all disinfectant residues.

• Additionally: Rinse all narrow and areas difficult to access, all cavities and lumen with a syringe several times (at least 3x) using each time 20 ml cold demineralized water.

• Dry product with a lint free wipe and sterile compressed air.

• Dry cavities, lumen and channels with sterile compressed air.

5.7. Automated Cleaning and Disinfection

Only use cleaning and disinfection devices (CDD) with proven efficiency according to EN ISO 15883.

Follow data of the manufacturer of the cleaning and disinfectant machine.

Only use agents suitable for medical devices made from metal and plastics with a pH-value between 5.5 and 12.3.

Recommendation: neodisher® mediclean forte (Dr. Weigert GmbH & Co. KG).

Apply program for thermal disinfection.

Follow instructions and data regarding program course and machine.

Product has to be stored safely and protected against mechanical damages during automated cleaning and disinfection.

Do not clean together with sharp edged or pointed objects.

Deposit product in a suitable rinsing basket.

Follow data for loading of the cleaning and disinfection device (CDD).

Use rinsing adapters for products with lumen and connect them according to the instructions in the User Manual of the cleaning and disinfection device (CDD).

Cleaning Program

Start the program course with following parameters:

- 1 Min. pre-rinsing with cold water
- Emptying
- 3 Min. pre-rinsing with cold water
- Emptying
- 5 Min. cleaning at 55°C with 0.5% alkaline cleaning agent
- Emptying
- 3 Min. neutralization with warm tap water (>40°C) and neutralizer (0.1% Neodisher® Z)
- Emptying
- 2 Min. interim rinsing with warm demineralized water (>40°C)
- Emptying

Disinfection Program

Automated thermal disinfection considering national requirements regarding AO value (see EN ISO 15883, AO value >3000).

• 5 Min. cleaning at 92°C +/-2°C

Drying

• 30 Min. at 90°C

• Remove rinsing adapter

At the end of program course, remove product and inspect it on remaining stain.

In case of residues, repeat automated cleaning and disinfection step as long as it takes until stain is no longer present.

Dry cavities and not sufficiently dried areas with sterile compressed air <2 bar.

Immediately after removal of product and immediately after additional drying on a clean place, put product in a single use sterilization packing (double packing) from paper or foil or put product in a sterilization container.

Respect requirements for sterilization packaging according to EN ISO 11607 and EN 868.

5.8. Sterilization

Only products that have been cleaned and disinfected are allowed to be sterilized.
Only apply steam sterilization in autoclave (fractioned pre-vacuum with sufficient product drying) for this product.

Adjust sterilization parameters:

• Minimum 134°C and maximum 137°C in saturated steam.

• Holding time at least 5 Min. until max. 20 Min.

• Drying in vacuum for at least 10 Min.

Sterilizer (Class B) according to valid national standards and regulation (e.g., EN 13060 or EN 285).

Recommendation: Sterilizer Class B, manufacturer: Tuttnauer.

Respect data of the sterilizer manufacturer regarding load, handling and drying times.

Exclusion:

Do not apply hot air, EO-gas, Radiation or Plasma for sterilization, or any other sterilization method for this product.

IMPORTANT:

Prior to usage, let product cool to room temperature.

It is the sole responsibility of the user to maintain the sterile condition of the product after the sterilization process.

In case the aforementioned chemicals and machines for cleaning, disinfection or sterilization are not available, it is the responsibility of the user to validate the procedure actually applied.

Also, if a sterilization method other than described above is applied, this deviating procedure has to be validated by the user accordingly.

5.9. Limitation of Reprocessing

The products are intended for single use only.

Therefore, reprocessing or reuse is strictly prohibited.

Only an undamaged product is allowed to be used.

6 Visual Inspection and Function Test

Prior to usage, the complete product, especially insulation and product tip have to be inspected for pressure points and damages.

A product exhibiting damages, pressure points or questionable condition is not allowed to be used and has to be replaced by a new one.

During application, tissue may adhere to the product, or sooting may be present on the distal end of the active electrode.

Such adhesions or sooting do not represent a reason for complaint and the product has to be exchanged by a new one.

Due to longer application time, mechanical forces or plasma seam or similar, the product may exhibit deformation or abrasion of the insulation material.

Also, such aspects do not represent a reason for complaint and the product has to be replaced by a new one.

Blockage of the suction channel (if applicable), does not represent a reason for complaint.

A product, that exhibits a blocked suction channel (if applicable), has to be replaced by a new one.

Prior to usage an electrical continuity test has to be done.

In case the product does not pass the electrical continuity test, it is not allowed to use the product any longer and has to be replaced by a new one.

7 Exclusion of repair and modification

Unauthorized repairs and modification (e.g. bending) are strictly prohibited.

Especially product, having a hook as active part, are never allowed to be bent.

This could lead to severe injuries of patient, user or third party.

Manufacturer's warranty immediately is terminated in case of modification or repair of the product.

8 Packaging, storage, transport

Store in a clean and dry environment, avoid direct sunlight.

Always handle the product with the utmost care when transporting, cleaning, disinfecting, maintaining, sterilizing and storing it.

This especially applies for fine tips and other sensitive areas.

Do not store or transport the product together with sharp edged or pointed objects.

Storage only in protective cases with individual compartments or individually sealed in foil.

Maintenance of the sterile condition after sterilization process is the sole responsibility of the user.

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REGER Medizintechnik GmbH, Gewerbestraße 10, 78667 Villingendorf, Germany
DR MEDICAL AG, Weissensteinstr. 28, Ziegelmatteareal Bau 6, 4503 Solothurn, Switzerland