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9 Hersteller

Dongguan QuanMed Equipment Mfg. Co., Ltd
Building 1, Dongyuan Industrial Center, Xiajiao Village
Caole Management Zone, Xiegang Town, Dongguan,
Guangdong – China
Tel: +86 769-82023123 / Fax: +86 769-82023163

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Sterile product for single use

HF-Pencil with Electrode Art.-No.: 8M50251500 (92649), 8M50251501 (92647).

Maximum rated voltage of accessory (U_{max}):

Description	U _{max}
Monopolar HF-Pencil with Electrode	5,0 kVp

See also label or catalog information.

In any combination with another electrosurgical accessory, the maximum rated voltage of the combination corresponds to the lowest rated voltage of the accessories used. (See also section "Intended Use")

See these instructions for use, the label, or the current product catalog for the maximum rated voltage of the product.

For use, read the entire IFU of this product and any used accessory as well as the HF generator and HF neutral electrode (monopolar application).

All requirements, safety notices, and warnings outlined in the respective IFUs must be strictly followed

This IFU is not intended for users in the USA.

2 Intended Use

Product intended for use in open or rather endoscopic surgery.

For cutting and coagulating of biological tissue.

Only to be used by skilled medical professional who have been introduced to the product, for example by this IFU.

Connect product to the monopolar output of the HF-generator by using the integrated cable of the device.

Activation is done by the buttons on the handle.

Footswitch can be used as well.

It is strictly prohibited to manipulate the electrode tip.

The product is supplied sterile and is for single use only.

It is strictly prohibited to reuse or reprocess the product.

Indication:

Product intended for use in open or rather endoscopic surgery.

Contraindication:

Product is not intended for direct contact with the heart or the central circulatory or central nervous system.

Application of high-frequency current may interfere with cardiac pacemakers and in vivo heart defibrillators, so affected patients must consult a cardiologist prior to the intervention.

It is not allowed the patient contacts any grounded metal objects during the intervention.

Connection and Activation:

Prior to application of monopolar products, it has to be safeguarded that an HF-neutral electrode has been applied correctly to the patient and that this device is connected properly to the respective HF-generator.

Activation is done using the buttons of the HF-handle or with the foot switch of the HF-generator.

Combination / Compatibility:

Prior to use, compatibility of the product with the intended HF-generator has to be verified.

In case of uncertainties, contact the manufacturer of this product or the manufacturer of the used HF-generator.

The frequency of the used HF-generator shall not exceed the maximum frequency of 4 MHz and not exceed the maximum rated voltage of the accessory (see section "Scope").

Follow instructions, safety notices and warnings included in the User Manuals or IFU of the used accessories and used HF-generator.

It is recommended to use a smoke evacuation system.

IMPORTANT:

Handle with utmost care.

This does not only apply for the duration of the surgery but also for the complete duration of storage and transport as well as during the process of connecting the product with the respective HF-handle, HF-cable and/or HF-generator.

This applies especially for the thin tips and other sensitive areas, for example the insulation.

Follow instructions and warnings in the IFU of HF-generator.

Side effects:

Electrosurgical interventions may lead to undesired damages or undesired burns of surrounding tissue.

Tissue damage due to extensive forces during butt dissection is possible.

It is the sole responsibility of the physician to include this aspect into their professional decision regarding the selection of treatment process.

All electrosurgical products potentially may lead to muscle stimulation during the application.

Improper use immediately will result in loss of warranty.

Liability for any damages incurred will not be accepted.

3 Safety Notices - WARNING

See these instructions for use, the label, or the current product catalog for the maximum rated voltage of the product.

Make sure the User Manual of the HF-Generator has been read carefully before using the product.

Prior to usage, a visual inspection and function test has to be done (see section "Visual Inspection and Function Test")

In case damages, deformation or similar is detected on the product, it is not allowed to use the device.

It has to be replaced by a new one.

At least one (1) sterile backup product has to be available for the application.

It is the responsibility of the user to determine the appropriate product size and product type according to their professional judgement and based on the patient's specific indication, preferred surgical technique and history, etc.

Prior to application, verify that the product is inserted firmly in the HF-handle/HF-cable. This must be done carefully, in order to avoid damages on the product and/or injuries of patient, surgical personnel or third party.

Excessive force can damage the product.

Therefore, the product has to be observed during the complete application.

Verify damage free sterile packaging and expiration date

Do not use a product that is included in a damaged sterile packaging or that is beyond the expiration date.

Maintenance of the sterile condition after sterilization process is the sole responsibility of the user.

Prior to usage a visual inspection has to be done for pressure points of damages on the product.

Damp inside connectors may cause flashover that could damage the connector or its insulation.

Prevent fluids from penetrating the connector during use.

Exclusion:

To plug in and unplug product, grasp it only on the plug part.

Do not activate the HF-generator during connection of the product.

Do not kink the product.

Do not use wound up or rolled up products.

As long as the product is activated, do not insert the tip in the surgical site and do not remove it when is activated.

Do not activate the product as long as it is in contact with metal objects and/or optics. During an electrosurgical intervention the patient must not come into contact with grounded metal objects such as for example surgical desk frames, instrument trays etc.

"Do not use in the product in the immediate vicinity of inflammable or combustible materials or substances. Active electrodes must not come into contact with inflammable or combustible substances such as disinfectants, gases etc. as otherwise there is a danger of fire or explosion due to heat or sparking."

Do not activate the product uninterruptedly over a longer period.

(Refer to section "General Safety Notices for HF-Technology")

After switching off the electrosurgical current, the product tip, still may be hot, so that this could lead to undesirable burns.

While the product is activated, do not insert or withdraw the tip from the surgical site.

Inadvertent activation of the product or moving the tip outside of the field-of-vision can lead to unintentional injury of tissue.

Do not use metal cannula, as these may damage the product insulation (capacitive coupling). This could lead to undesirable burns.

It is recommended to use cannula systems made of plastics.

It is the responsibility of the user to apply low performance settings on the HF-Generator in order to achieve the envisaged effect for the respective intervention.

During the surgery, mechanical forces may lead to deformation or signs of wear of the product.

If abrasion of the product is detected during the application, the product has to be exchanged by a new one.

Abrasion, adhesion of tissue, discoloration, soothing etc. do not represent a reason for complaint and do not permit claiming the manufacturer's warranty.

4 General Safety Notices for HF-Technology (Excerpt)

In addition to the acknowledged benefits of HF-Technology, the application includes several risks that have to be attended:

Improper use and non-observance of the IFU can lead to unintentional burns of the patient, user or third party.

Continued further education of the surgical personnel is recommended.

a) Environment

Pay attention that no flammable substances (anesthetics, oxidizing gases, endogenous gases etc.) are present in immediate vicinity as otherwise there is danger of explosion.

Only use non-flammable cleaning or disinfection agents.

All oxygen connections must be tight and leak-proof during the procedure.

b) Patient Positioning and Preparation

Ensure proper patient positioning, i.e., only use insulating surgical drapes that are dry, absorbent and liquid-tight.

Isolate conductive surfaces and contact points towards the patient.

Dry tissues of pulp are necessary to use for skin and breast folds as well as between the extremities.

Prior to application, remove any liquids that potentially accumulated in body cavities

Only use non-flammable disinfectants.

Do not use alcohol-based lotions.

Only use non-conductive rinsing fluid, if medicinally possible.

Attend the requirements on irrigation fluid for monopolar and bipolar products.

Remove all body jewelry from the patient.

Putting a band aid over the body jewelry is not sufficient

c) Connections

Prior to application, make sure the used HF-handle/HF-cable is connected properly to the HF-generator and make sure the correct power and performance setting is adjusted and displayed.

It is the responsibility of the user to apply low performance settings on the HF-Generator in order to achieve the envisaged effect for the respective intervention.

d) HF-Neutral Electrode for Monopolar Application

In case of monopolar application, select a HF-neutral electrode suitable for the patient, apply it correctly and connect it properly with the respective HF-generator.

Follow all instructions for proper application of the HF-neutral electrode, incl. patient protection and patient monitoring of the HF-neutral electrode and all further provisions, safety notices and warnings included in the IFU of the HF-neutral electrode.

e) Patient Reactions

All electrosurgical devices potentially can cause muscle stimulation during the application.

The design of this product minimizes the risk of this undesirable effect.

Nevertheless, muscle stimulation can lead to an unexpected movement of the patient in the surgical field.

f) Handling HF-Accessory

Make sure the accessory is compatible.

Do not touch the instrument tip during the complete application.

As long as the product is not applied, place it on a dry, clean and non-conductive and well-visible surface, that is not in contact with the patient.

Never store product on the patient.

Unintended activation of the product can lead to unintentional burns or other injuries of patient, user or third party.

Never wrap cables around the patient and never lay cables over the patient.

Apply only short activation times.

Keep on longer breaks between activation phases.

Only adjust low power settings.

g) Completeness of the System

At the end of the surgery, confirm completeness of the system.

5 Visual Inspection and Function Test

Prior to usage, verify the sterile packaging is free of damage and the expiration date is valid.

Do not use a product that is included in a damaged sterile packaging or that is beyond the expiration date.

It has to be replaced by a new one.

Never lay this product on the patient or in their direct vicinity.

Prior to application, verify the correct function of the buttons:

Blue button: COAG / Yellow button: CUT

If buttons do not work, use footswitch.

Note: If the cable conductor is broken or the electrical continuity of the cable is interrupted in any other manner it can cause arcing in the patient return system or in the active circuit.

This can lead to the patient suffering burns or the outbreak of a fire.

If damaged, do not use the product.

Prior to each usage, the complete product, especially insulation and product tip have to be inspected for pressure points and damages.

A product exhibiting damages, pressure points or questionable condition is not allowed to be used and has to be replaced by a new one.

During and after application, tissue may adhere to the product, or soothing may be present on the distal end of the active electrode.

Such adhesions or soothing do not represent a reason for complaint and the product has to be exchanged by a new one.

Due to longer application time, mechanical forces or plasma seam or similar, the product may exhibit deformation or abrasion of the insulation material.

Also, such aspects do not represent a reason for complaint and the product has to be replaced by a new one.

Blockage of the suction channel (if applicable), does not represent a reason for complaint.

A product, that exhibits a blocked suction channel (if applicable), has to be replaced by a new one.

Prior to usage an electrical continuity test has to be done.

In case the product does not pass the electrical continuity test, it is not allowed to use the product any longer and it has to be replaced by a new one.

6 Exclusion of Repair and Modification

Unauthorized repairs are strictly prohibited.

Unauthorized modification is strictly prohibited and lead to an immediate loss of manufacturer's warranty.

It is not allowed to bend product or to manipulate the product in any other manner.

Especially product having a hook as active part is never allowed to be bent.

This could lead to severe injuries of patient, user or third party.

The product is delivered sterile and for single use only.

It is not allowed to reuse or reprocess the product.

See also labelling.

7 Packaging and Labelling

Only use product in closed original packaging with an undamaged labelling.

In case of uncertainties, contact the manufacturer.

Only use product of an undamaged sterile packaging and valid expiration date.

Do not use product, that exhibits a damaged sterile packaging or that is beyond the expiration date.

See also labelling.

8 Storage

Store product in a clean and dry environment protected from direct sunlight.

Store product in the closed original packaging.

Maintenance of the sterile condition after sterilization process is the sole responsibility of the user.

Do not use the product any longer if the expiration date is expired or that has a damaged sterile packaging.

9 Manufacturer

Dongguan QueenMed Equipment Mfg. Co., Ltd

Building 1, Dongyuan Industrial Center, Xiajiao Village

Caole Management Zone, Xiegang Town, Dongguan,

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info@dr-medical.ch | www.dr-medical.ch

13 Warranty

Liability and warranty are excluded for all products that have been modified in any way, not applied according to their intended use, or that have been handled or applied improperly, or in case of any other deviation from the instructions set forth in this IFU.

Furthermore, the manufacturer denies any liability for any accidental, intentional damage or for a damage or loss arising out of handling or application of the product.

Additionally, all liability and warranty are omitted in case our product was repaired by a company that has not been authorized by us.

Unauthorized repairs are strictly prohibited.

14 Return

Returned products only are accepted, if they are marked as "hygienically safe" or "not contaminated" and have been packed safely for the transport.

Use our return sheet for any return of products.

15 Disposal

Disposal of the product, their packaging material as well as any accessory has to be done according to the applicable country specific requirements, regulations and laws.

Additionally, the applicable requirements of the respective clinical place(s) in regards to disposal of medical devices have to be followed as well.

16 Regulatory Remark

Due to regulatory requirements, we would like to inform you that any serious incidents related to this medical device must be reported to the manufacturer as well as to the competent authority of the EU Member State where the user and/or patient resides.

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18 About these Instructions for Use

The instructions for use must be kept in a location that is easily accessible to all users throughout the entire period of the product's use.

For a current revision of these IFU, please contact medimex GmbH or visit the medimex Homepage for a download of this document:

www.medimex.de/downloads

Subject to change.

This IFU is not intended for users in the USA.

©2024. REGER Medizintechnik GmbH All rights reserved. All company names, service marks, trademarks, product names and logos appearing on this document are the property of their respective owners and shall be treated accordingly.