

The ^{new} mediPlugs⁺



New Design – Improved Quality!

Our feet carry us 130,000 kilometers through life, enough to circle the globe nearly three times. They bear a load of 2,500 tons every day – for an average of 78 years. With 26 bones, 30 joints, 60 muscles, 100 ligaments, and 200 tendons, they are a true masterpiece of the human body. Are you impressed? So are we! medimex stands for decades of experience in developing the highest-quality surgical shoes.



The orthopedic insole provides support for your feet.



The special ventilation system keeps your feet fresh.



The excellent tread pattern ensures a secure footing.



The snap-fit system connects the sole and the shoe during refurbishment.

mediPlogs⁺: The health benefit

Always in balance thanks to orthopedic insoles

- the ergonomic shape of the insole and its textile, antibacterial surface provide exceptional support
- the insole supports the transverse arch, prevents the foot from sliding forward, and stimulates the toes' gripping function
- the contoured edge of the insole centers the foot and keeps it in the correct position

Air-conditioned via a special ventilation system

- keeps your feet comfortable
- the design of the side ventilation holes largely prevents liquids from entering the shoe

The extra touch of safety & service

- the special design of the heel, toe, and rounded edge of the outsole reduces the risk of tripping and stabilizes your gait
- the snap-in system secures the insole inside the shoe and allows it to be attached to the outside of the shoe during processing – no searching, no sorting

Machine washable up to 60°C

- mediPlogs⁺ are machine washable up to 60°C and can be cleaned with any machine-safe disinfectant

Made of antistatic PU

- CFC-free foam
- ensures safety and protects the environment

Meets the applicable essential health and safety requirements

- complies with the requirements of the harmonized standard EN ISO 20347:2022 and EN ISO 20347:2022/A1:2024
- PPE protection class: OB A E
- PPE risk category: II

Customization service for mediPlogs⁺ with heel straps

- the mediPlogs⁺ with heel straps can be customized with your name or department code upon request

Name	Colour	Size	Pair/ packing unit	Item no. without heel strap*	Item no. with heel strap*	Item no.*
mediPlogs ⁺ OP-shoe	green	36 – 47	1	8M1374 __	8M1384 __	
mediPlogs ⁺ OP-shoe	blue	36 – 47	1	8M1375 __	8M1385 __	
mediPlogs ⁺ OP-shoe	fuchsia	36 – 47	1	8M1376 __	8M1386 __	
mediPlogs ⁺ OP-shoe	light blue	36 – 47	1	8M1377 __	8M1387 __	
Replacement insole	light blue	36 – 47	1			8M0345 __

*Please include the size with the item number when ordering.

